




[www.ajmostar.com](http://www.ajmostar.com)

## Orthodontic Instruments Catalogue



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## WEINGART UTILITY PLIERS



30-502T

WEINGART UTILITY PLIER



60-501T

WEINGART UTILITY PLIER



60-501

WEINGART UTILITY PLIER



20

WEINGART UTILITY PLIER

AJMO STAR  
SURGICO



**Ajmo Star Surgico** Showall(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments

## UTILITY PLIERS

### HOW PLIERS



**60-306T**

HOW PLIER  
Tips are 1/8" (3.2mm) diameter.



**60-306**

HOW PLIER  
Tips are 1/8" (3.2mm) diameter.



**60-307**

HOW PLIER-SMALL  
Tips are 3/32" (2.4mm) diameter.



**60-308**

HOW PLIER-OFFSET  
Tips are 1/8" (3.2mm) diameter.



**60-219**

LIGATURE UTILITY PLIER



The width between two handles measures about 1cm less than the regular How plier and it allows rotating of the plier easily.



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments

## WEINGART UTILITY PLIERS



30-502T

WEINGART UTILITY PLIER



60-501T

WEINGART UTILITY PLIER



60-501

WEINGART UTILITY PLIER



20

WEINGART UTILITY PLIER



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## BENDING PLIERS

### BIRD BEAK PLIERS



**30-314T**

BIRD BEAK PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .032"



**60-314**

BIRD BEAK PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .030"

### LIGHT WIRE PLIERS



**30-315C**

LIGHT WIRE PLIER WITH CUTTER  
-GROOVED SQUARE TIP  
MAXIMUM BENDING & CUTTING CAPACITY:  
Wires up to .020"



**30-315T**

LIGHT WIRE PLIER  
-GROOVED SQUARE TIP  
MAXIMUM BENDING CAPACITY:  
Wires up to .020"



**60-315**

LIGHT WIRE PLIER  
-GROOVED SQUARE TIP  
MAXIMUM BENDING CAPACITY:  
Wires up to .018"



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Slalkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments



**60-216**

ELASTIC PLACING PLIER-HOOK TIP



**60-215R**

LIGATING PLIER WITH ROUND HEAD



With a round head, you can hold a ligature wire certainly.  
A time for ligating should be shortened. The cross spring  
may reduce fatigue with very smooth opening and closing  
of the instrument.



**60-217**

LIGATING PLIER



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## BRACKET PLACERS



**60-815**

BRACKET PLACING TWEEZERS



**60-890**

BRACKET SETTING TWEEZERS



DIAMOND TWEEZERS

The diamonds are electrodeposited inside the point of the tweezers. Being not to slide, ligature wires or brackets can be carried into a mouth properly.

## BRACKET POSITIONING GAUGE

**60-830**

BOON BRACKET POSITIONING GAUGE  
positioning height : 3.5mm/4.0mm/4.5mm/5.0mm



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## BANDING PLIERS

### BAND REMOVING PLIERS



**60-983**  
REPLACEMENT PEEK TIP(3pk.)

**60-104T**  
POSTERIOR BAND REMOVING PLIER



**60-983**  
REPLACEMENT PEEK TIP(3pk.)

**60-104TL**  
POSTERIOR BAND REMOVING PLIER-LONG



**60-985**  
REPLACEMENT ALUMINUM TIP(3pk.)

**60-104**  
POSTERIOR BAND REMOVING PLIER



**966**  
REPLACEMENT ALUMINUM TIP(3pk.)

**30**  
POSTERIOR BAND REMOVING PLIER-LONG



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)





# orthodontic instruments

## ARCH FORMING INSTRUMENTS

### ARCH FORMING PLIERS



**60-303T**

TWEED ARCH BENDING PLIER

Beaks are .060" wide, parallel at .020" opening

MAXIMUM BENDING CAPACITY:

■ Rectangular Wires up to .022" x .028"



**60-303**

TWEED ARCH BENDING PLIER

Beaks are .043" wide, parallel at .020" opening

MAXIMUM BENDING CAPACITY:

■ Rectangular Wires up to .022" x .028"



**60-006**

HOLLOW CHOP PLIER

MAXIMUM BENDING CAPACITY:  
Wires up to .022"



**60-223TG**

ARCH FORMING PLIER-GROOVED

Grooves of .016", .018" and .022"



**60-223T**

ARCH FORMING PLIER-NON GROOVED



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## DEBONDING PLIERS

### BRACKET REMOVING PLIERS



PRORTHO



**P-200**

**BRACKET REMOVING PLIER**  
The blades can catch the bracket's edge by sliding on the surface of the tooth and come beneath the attachment without efforts. The blade edge comes to touch at point, not a line contact, it reduces resistance for easy control.



**60-334T**  
DIRECT BOND ATTACHMENT REMOVING PLIER-ANGLED



**60-335T**  
DIRECT BOND BRACKET REMOVING PLIER-STRAIGHT



**TS-25 (TASK SPECIAL)**  
DIRECT BOND BRACKET REMOVING PLIER-STRAIGHT



**TS-26 (TASK SPECIAL)**  
DIRECT BOND ATTACHMENT REMOVING PLIER-ANGLED



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments

## CRIMPING PLIERS

### HOOK CRIMPING PLIERS



**PRORTHO**



Pressing mechanism  
(front view)



Installing the hook

**P-220**

#### HOOK CRIMPING PLIER

The beak that has unique grip form pressures a hook firmly on the wire to lock securely. The flat handle of PRORTHO is/ears slots on the palm.



**60-446**


HOOK CRIMPING PLIER  
For quick and precise installation of hook on the archwire.

### TUBE CRIMPING PLIER



**60-319**

TUBE CRIMPING PLIER



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments



**60-817**  
LIGATING & DIRECTOR



**60-818**  
LIGATURE DIRECTOR-DOUBLE END



**60-819**  
LIGATURE DIRECTOR-SINGLE END



**60-816**  
LIGATURE DIRECTOR



**SH DIAMOND LIGATURE DIRECTOR**  
The diamonds are electrodeposited in the point. Being not to slide, Wire can be placed properly.





# orthodontic instruments

## LIGATING INSTRUMENTS

### LIGATING INSTRUMENTS



**60-219**

LIGATURE UTILITY PLIER

The newest concept in utility pliers for ligating operation with module or ligature as well as the other multiple services How pliers offer.



Comparing with the standard How Plier, this new model forms a narrow oval between the closed handles. The width between two handles measures about 1cm less than the regular How plier and it allows rotating of the plier easily.



operation with module



operation with ligature wire



**60-218**

ELASTIC PLACING PLIER



**N-801**

MOSQUITO FORCEPS-STRAIGHT



**N-802**

MOSQUITO FORCEPS-CURVED



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)



# orthodontic instruments

## BENDING PLIERS

### LOOP FORMING PLIERS



**60-074**

YOUNG PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .040"



**60-213**

JARABAK PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .020"



**60-305**

LOOP FORMING PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .022"



**60-304**

TWEED LOOP FORMING PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .022"



**60-304S**

OMEGA LOOP FORMING PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .022"



**60-304TR**

TWEED LOOP FORMING PLIER  
MAXIMUM BENDING CAPACITY:  
■ Rectangular Wires up to .022" x .028"



**60-990**

Replacement Tip (1/pk.)  
(Two Rounds Size: .043" .063")

**60-971**

HEX WRENCH (1/pk.)



**Ajmo Star Surgico**

Showali (East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments



**P-500**  
LINGUAL BRACKET REMOVING PLIER



## BOND REMOVING PLIERS

PRORTHO



**P-210**  
BOND REMOVING PLIER



PRORTHO




**60-983**  
REPLACEMENT PEG TIP (3/pk.)



**P-211**  
REPLACEMENT METAL TIP (1/pk.)  
both ends can be used



**60-971**  
HEX WRENCH (1/pk.)



**60-332T**  
BOND REMOVING PLIER



**60-9871**  
REPLACEMENT METAL TIP (1/pk.)



**60-9872**  
REPLACEMENT TUBE (15mm 1/pk.)



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments

## LIGATING INSTRUMENTS

### LIGATING INSTRUMENTS



**60-101C**  
LIGATURE TYING PLIER  
-COON TYPE (SAFETY TIPS)



**60-102**  
LIGATURE TYING PLIER



**60-225T**  
LIGATURE FORMING PLIER



**N-215**  
MATHIEU TYPE







# orthodontic instruments



**60-224T**

STOP PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .025"

CINCH BACK PLIER



**60-900**

TIE-BACK PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .025" (Ni-Ti Only)



QUAD HELIX BENDING PLIERS



**60-202T**

QUAD HELIX BENDING PLIER  
MAXIMUM BENDING CAPACITY:  
Wires up to .036" (0.9mm)

**60-970**  
Replacement Tip (1/pk.)  
Round Size: .063" (1.6mm)

**60-971**  
HEX WRENCH (1/pk.)  
©Details the list of page 139 viewing



**60-201L**

ANGLED THREE JAW PLIER-LEFT  
MAXIMUM BENDING CAPACITY:  
Wires up to .036" (0.9mm)  
Used for a small adjustment of quad helix  
or lingual arch. Narrow beaks allow an easy  
access to lower jaw.



**60-201R**

ANGLED THREE JAW PLIER-RIGHT  
MAXIMUM BENDING CAPACITY:  
Wires up to .036" (0.9mm)  
Used for a small adjustment of quad helix  
or lingual arch. Narrow beaks allow an easy  
access to lower jaw.





# orthodontic instruments



60-430CK



60-302

For bending of multiple edgewise arch wires.

MAXIMUM BENDING CAPACITY:

■ Rectangular Wires up to .016" x .022"

Beaks are parallel at .020" opening

MAXIMUM BENDING CAPACITY:

Wires up to .022"



60-303T

Beaks are .060" wide, parallel at .020" opening

MAXIMUM BENDING CAPACITY:

■ Rectangular Wires up to .022" x .028"



60-303

Beaks are .043" wide, parallel at .020" opening

MAXIMUM BENDING CAPACITY:

■ Rectangular Wires up to .022" x .028"



**Ajmo Star Surgico**

Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: info@ajmostar.com, sales@ajmostar.com



# orthodontic instruments

## BAND CONTOURING PLIERS



**60-208**  
BAND CONTOURING PLIER



**60-212**  
CROWN CONTOURING PLIER



**60-317**  
CROWN CRIMPING PLIER



**60-317S**  
CROWN CRIMPING PLIER-SLIM



**Ajmo Star Surgico** Showali(East.) P.O. Cantt Sialkot-51310, Pakistan

Tel: +92 52 4273146 Fax: +92 52 4298566 Email: [info@ajmostar.com](mailto:info@ajmostar.com), [sales@ajmostar.com](mailto:sales@ajmostar.com)